



www.spermprocessor.com



Sperm Processor Pvt. Ltd.

6, Welcome Nagar, Garkheda,
Aurangabad (MS) - 431005, India.

Ph.: +91 240 6603800

Fax: +91 240 2341694

Email: info@spermprocessor.com

α -Glycosidase

bio-marker of epididymis activity

REF] SP/SFT/ α -011



IVD



Turnaround time for test: 140 min



Store at: 2°C - 8°C after receiving

Index

Concept	01
Specimen Preparation	02
Special Instructions	02
Kit Contents	03
Equipments	04
Disposable Materials	04
Bio-chemistry Analyser Settings	05
Procedure	06
Preparation of Serial Standard Graph ...	07
Examination	08
Result	10
Precautions.....	11
Safety & Environment	11

α - Glycosidase

Bio-marker of epididymis activity

CONCEPT

Poor quality semen may result from testicular production of abnormal sperm or from post testicular damage of sperm in the epididymis or the abnormal ejaculate from accessory gland. Male accessory sex organs made of epithelial / mesenchymal components, require 'Androgen' for proliferation & maintenance of their functions.

The epididymis makes a significant contribution to the development of a fertile sperm, through four classical functions, sperm maturation, sperm transport, sperm concentration & sperm storage.

α -Glycosidase represents as Bio-marker of Epididymis as its major contribution is from Epididymis. The evaluation of α -Glycosidase is regarded as a reliable parameter of Epididymal potency.

Post-testicular maturation of sperm taking place in the epididymis is now well recorded, makes a significant contribution to the development of fertile ejaculate.

Differential diagnosis of Azoospermia only the determination of the α -Glycosidase is useful.

The reduction of α -Glycosidase in infertile male is to be attributed to reduced sperm count.

1 Specimen Preparation

- Semen sample is collected with :
 - **Abstinence period** of **2-7days**.
 - **Ideal collection** through **masturbation** in sterile container.
 - **Non-spermicidal polyurethane semen collection pouch (Sperm Collect™)** can be used when required.
- Semen sample is allowed to liquefy and then well mixed for performing test.
- **Sperm-free Seminal Plasma Preparation :**
 - **Note down** semen volume (upto one decimal).
 - **Centrifugation** (with Androspin™) of semen sample (liquefied, well mixed) at **3000 rpm** for **10 - 20 min**.
 - **Aspirate supernatant** to obtain seminal plasma & **leave the pellet**.
 - Can be stored at **-20°C** for **12 months**.


Special Instructions :

- **Hyperviscous** semen sample should be processed to bring towards normal viscosity. (**Viscosity-CH™** or **Viscosity-BR™** kit can be used)
- **Frozen semen** plasma must be thawed at 37°C (with Sperm Warmer™) before performing test.

2 Kit Contents

- Reagent 1 (R1) Reaction buffer : 10 mL
- Reagent 2 (R2) Substrate solution (Conc.) : 0.4 mL
- Reagent 3 (R3) Inhibitor solution : 10 mL
- Reagent 4 (R4) Stopping buffer : 120 mL
- Reagent 5 (R5) Standard dilution buffer : 60 mL
- Reagent 6 (R6) PNP Stock Standard (05 mM): 02 mL

Kit Content Layout Diagram :



Storage Conditions :

- The kit should be stored in dark at 2°C - 8°C after receiving.
- Bring all the reagents to room temperature before use.
- Once opened, store reagents in the fridge protected from light.
- Expiry date is printed on the out side of the box.

REQUIRED BUT NOT PROVIDED IN KIT

- Controlled Temperature 37°C Dry bath (Sperm Warmer™ / Water bath)
- Centrifuge Machine (Androspin™)
- Bio-chemistry Analyser (Androchem Analyser™)
- Pipettes Set
- Stopwatch
- Microtip Box
- Test Tube Stand

REQUIRED BUT NOT PROVIDED IN KIT

- Hand gloves
- Semen Collection Container
- Non-spermicidal Semen Collection Pouch (Sperm Collect™)
- Microtips
- Pasteur Pipettes
- Test Tubes
- Micro Tubes

Assay Parameters For Programming :

- | | |
|-----------------------------|------------------------------------|
| 01. Test Name | : Alpha Glycosidase |
| 02. Mode | : Differential with multi standard |
| 03. Primary wavelength | : 405 nm |
| 04. Temperature | : 37° C |
| 05. Aspiration Volume | : 500 µL |
| 06. Lag time | : 05 Sec |
| 07. Blank | : No |
| 08. Q. C. | : No |
| 09. Standard | : 05 |
| 10. Concentration | : 0 , 10, 50, 100 & 150 µM |
| 11. Norm | : Yes |
| 12. Reagent linearity | : 250 mIU/ mL |
| 13. Read sec | : -- |
| 14. K factor | : 0.484 |
| 15. Unit | : mIU / mL |
| 16. Reaction Time | : 120 min |
| 17. Reagent Blank Abs. max. | : 1.5 |

AndroChem[™]
Analyser



(Mfg. by Sperm Processor Pvt. Ltd.,
Aurangabad, India)

- Is preprogrammed with above mentioned settings
- No additional programming is required.

- Step 1 :**
- **Note down semen volume** (upto one decimal).
 - **Bring** all the reagents to room temperature before use.

- Step 2 :** Take **two** micro centrifuge tubes & label them as '**Blank**' (**B**) & '**Test**' (**T**).
(With appropriate **Lab ID No.**)

- Step 3 :** Take the reagents & put in appropriate tubes with help of air displacement pipettes as per the following table.



Reagents	Test Tube	
	Blank (B)	Test (T)
R1 – Reaction buffer	140 μ L	----
R3 – Inhibitor solution	----	140 μ L
R2 – Substrate solution	03 μ L	03 μ L
Seminal Plasma	20 μ L	20 μ L

- **Incubate tubes at 37°C (Dry Bath) for 2 hrs.**
- **After 2 hrs add following reagent (R4)**



R4 - Stopping Buffer	1000 μ L	1000 μ L
----------------------	--------------	--------------

NOTE :

- Each test preparation require own blank.
- Avoid common blank, whenever there are more than one test.

- Step A :**
- Take a test tube & prepare **Solution - A (Diluted PNP Std.)** with the help of following as per following table.

R6 (PNP Stock Std.)	80 μ L
R5 (Std. Dilution Buffer)	1920 μ L
Solution - A (Diluted PNP Std.)	2000 μ L


- Mix well & use solution for **Step B**

- Step B :**
- Take 5 new test tube & label them as **Std.1 (S 1), Std. 2 (S 2), Std. 3 (S 3), Std. 4 (S 4) & Std. 5 (S 5)**
 - Add reagents to respective test tubes as per table given below & mix well.

Test Tube	Step A Solution (Diluted PNP Std.)	Reagent 5 (R5) (Std. dilution buffer)
Std. 1 (0 μ M)	----	1000 μ L
Std. 2 (10 μ M)	050 μ L	0950 μ L
Std. 3 (50 μ M)	250 μ L	0750 μ L
Std. 4 (100 μ M)	500 μ L	0500 μ L
Std. 5 (150 μ M)	750 μ L	0250 μ L

- Step C :**
- First read absorbance of **Std. 1 (S 1)** by using **distilled water** as **blank** .
 - Read serially absorbance of **Std. 2 (S 2), Std. 3 (S 3), Std. 4 (S 4) & Std. 5 (S 5).**

Step D : Plot **Serial Standard (Std) Graph** as shown below by plotting concentration of Std. on X axis & absorbance on Y axis.



- If **Androchem Analyser** is used, results are displayed automatically. If required take print out.

NOTE : Save this Std. Graph & use for test results.

8 Examination

- Read the absorbance by using Bio-chemistry Analyser (at 405 nm), within 60 minutes.
- Read the absorbance of Blank.
- Prepare the serial standard graph.
- Read the absorbance of test.

NOTE : For each test prepare new blank & avoid using common blank.

Calculations

α-Glycosidase in semen per ejaculate

1. α-Glycosidase in semen (mIU/mL)

A. If Biochemistry Analyser is used :

- α-Glycosidase in semen (**X**) mIU / mL :

$$[\text{Abs. T} - \text{Abs. B}] \times (\text{Std. curve factor}) \times 0.484$$

$$= \mathbf{X \text{ mIU / mL}}$$

Note : Correction factor 0.484 is calculated on the basis of sample dilution & incubation time.

- Final Result
 α-Glycosidase in semen sample = **X** mIU/mL

B. If AndroChem Analyser is used :

*AndroChem*TM
Analyser



(Mfg. by Sperm Processor Pvt. Ltd.,
Aurangabad, India)

- Is preprogrammed with above mentioned settings.
- No additional programming is required.
- Results are displayed as **X** mIU / mL

Quantitative Estimation of α -Glycosidase in Seminal Plasma

(Biomarker - Epididymis)

Volume : ____ mL

Result : ____ mU / mL

: ____ mU / ejaculate

Normal Reference Range :

≥ 20 mU / ejaculate

(As per fifth edition of WHO laboratory manual for examination and processing of human semen).

Limitations :

- This test provides presumptive quantitative information of α -Glycosidase in seminal plasma
- This parameter should be analyzed by a specialist.
- The result should be evaluated taking into account all clinical & laboratory findings related to the same sample.

- All patient samples & reagents should be treated as potentially infectious & the user must wear protective gloves, eye protection & laboratory coats when performing the test.
- The kit should be discarded in a proper biohazard container after testing.
- Do not eat, drink or smoke in the area where specimens & kit reagents are handled.
- Do not use beyond the expiration date which appears on the package label.
- It is recommended to use of gloves & face mask.

- Do not release the products used into the environment. Follow centre guidelines for the storage & disposal of toxic substances.
- Biological samples must be handled as potentially infectious.

Description of Symbols



consult instructions of use



product reference



lot number



use by



manufacturer



health surveillance device
for in-vitro diagnostic



contains sufficient for 'n' tests



temperature limitation



keep dry



CE mark (Conformité Européene)

Accreditations & Registered Certificates

- **ISO 13485 : 2003** Certified
- **CE** Certified
- **GMDN** Registered
- **US FDA** Registered

For more information & procedure videos



www.spermprocessor.com/sft-alpha-glycosidase.html



www.youtube.com/watch?v=2dfXIVBjC-4

all logos are trademarks
and sole proprietary of Sperm Processor Pvt. Ltd.