





Sperm Processor, an ISO 13485:2003 company, presents **Sperm 360°** - an **end to end solution** for andrology labs.

All our products are **CE accredited**.

Each product is thoughtfully designed and precisely made to increase performance quotient of your Lab.

index

Semen Analysis & SFT Kits	1
Instruments & Accessories	19
CASA Setup	29
Plasticware & Disposables	35
Lab Setup Solutions	39



360° Solutions

that redefine the way to serve Lab Andrology...

widest product gamut

turnkey setup solutions

franchisee model for lab services

index

Semen Analysis & SFT Kits	
Instruments & Accessories	1
CASA Setup	2
	ľ
Plasticware & Disposables	1
Lab Setup Solutions	3

2

3

Routine Semen Analysis

Viscosity

a reagent for hyper viscous semen



needs different capability.

now you can check each of it

with our SFT Kits

from production to sperm-ovum interaction sperm gets exposed to different environments; for each environment it



Viscosity-CH: Chemotrypsin Solution Viscosity-BR: Bromolene Solution Product Code SP/SFT/VS-001-A SP/SFT/VS-001-B



Semen Analysis & SFT Kits

ready to use kits

 $supported\ with\ specially\ designed\ accessories$

every kit come with video of test procedures

all products are

(

accros

Vitality

evaluation of sperm membrane structural integrity



HOS

evaluation of sperm membrane Unctional & structural integrity



Product Code SP/SFT/V-003-A

Vitality I - Eosin -Nigrosin Solutions

Product Code SP/SFT/H-004

Semen pH reagent strips for determination of pH of semen



Leukocytes

detection of leukocytes by using peroxidase test





Product Code SP/SFT/PH-002

Product Code SP/SFT/L-006



Morphology-D

evaluation of sperm morphology by Diff-Quick Staining method



Product Code SP/SFT/MD-005-A

Morphology-P

evaluation of sperm morphology by Pap Staining method

introducing spray method



Product Code SP/SFT/MP-005-B

10

1



in seminal plasma

biomarker of epididymis activity



also see

Androchem Analyser Pg 24

quantitative estimation of Fructose in seminal plasma





biomarker of iinal vesicle activity

Product Code SP/SFT/F-012

Product Code SP/SFT/α-011

Sperm Content Test



quantitative estimation of

Zinc

in seminal plasma

biomarker of orostate activity



Nuclear Protein Assay

to assess nuclear maturity



Product Code SP/SFT/NP-008

Product Code SP/SFT/Z-013



DNA Fragmentation

sperm chromatin dispersion test



Product Code SP/SFT/DNA-009



DNA Fragmentation kit includes hydrophilic agarose coated slide & cover glass

also see

we have specially developed accessories to give you complete DNA Fragmentation Solution



DNA warmer - Pg 23



DNA slide immersion tray - Pg 24



Antisperm Antibody - IgG

MAR Test (mixed antiglobulin reaction)



fresh semen can be use for testing

Antisperm Antibody IgG - Latex Particles coated with IgG and Antihuman IgG Antiserum

Product Code SP/SFT/ASA-007-A

Antisperm Antibody - IgA

MAR Test (mixed antiglobulin reaction)



resh semen can be used for testing

Antisperm Antibody IgA - Latex Particles coated with Antihuman IgA Antibody

Product Code SP/SFT/ASA-007-B





Instruments & Accessories

Practical Innovations for Practical Andrology

 $ergonomic\ designs\ |\ greater\ operational\ functionality\ |\ improves\ test\ results$

all products are



credited

Sperm meter

SpermMeter analyses **sperm concentration** and with respect to it, quantitative & qualitative assessment of **sperm motility**.

importantly, you get all these results in a same focus of microscope

no need of semen dilution

 $10 \mu \; depth \; avoids \\ overlapping \; of \; sperms$

defined micrometric grid for easy calculations



Product Code SP/INT/001-A

Sperm

SpermWarmer is a microprocessor based **temperature control system** to maintain

specimens at 37° C temperature



provision to hold all the plasticware, & disposable required for the test:

Product Code SP/INT/004-A

Slide Warmer

Slide Warmer is a microprocessor based temperature control system to heat-dry prepared smears & stained smears at 40°C-50°C



ensures uniform heat transfer

saves time

Product Code SP/INT/004-B



DNA Warmer is a dual temperature control system designed to expose the agarose tubes at 70°C or at 37°C in separate blocks of instrument



displays set & current temperature

compatible for 0.5 ml conical bottom agarose tubes

Product Code SP/INT/004-C

DN Slide Immersion Tray

specially designed tray for DNA Fragmentation Test to immerse the agarose coated slides completely with **minimal quantity of reagents**

saves reagent

reduces cost per test

Product Code SP/ACC/003-B

Glass Slide Staining Tray

a tray to **maintain** slide surface in exact **horizontal position**for layering & maintaining the reagents on it
during staining procedure

only the required portion of reagent is utilise

no chances of change in pH of the rest of reagent

improves staining quality

Product Code SP/ACC/003-A

AndroChem Analyser

Androchem Analyser is pre-programmed to perform **quantitative estimation of** the **biomarkers** like Fructose, Zinc, α - Glycosidase



instant report prin

R-32 port to transfer the results to computer

facility to program new test

Product Code SP/INT/008



a **centrifuge machine** designed to **work at 37°C**, that is same as the in-vivo temperature, to maintain the natural functionality of sperm



extremely sturdy

Product Code SP/INT/006-A

Sperm visičn

phase contrast microscope -

10X & 20 X phase optics with 40X & 100X bright field optics with all standard configurations





Computer Assisted Semen Analyser

Product Code SP/INT/005-A

Sperm secto CASA & You

for the first time ever a CASA that includes You...

introducing Expert Mode

where the technology helps you go beyond technology



Product Code SP/INT/002

patent pending technology

SpermSoft is par image processing It's is not just a software but a complete stand-alone system that will surely be a backbone of your andrology lab...

SpermSoft is developed with a bunch of exciting features so that you can work with ease & convenience and be always confident about the results





traceability with har code



image & video reports



ease of touch scree

stay connected through wi-fi

SpermSoft comes with wi-fi connected tablet, so now you can monitor or give your expert input to lab technicians right from your desk



3

world's most comprehensive CASA

E	xpert Mode	Auto Mode	
Routine Semen Analys	is 🗸	/	
Morphology	/	/	
Vitality	/	✓	the stand-alone system comes with
HOS	~	NA	negative phase contrast
Nuclear Protein Assay	/	✓	microscope & all accessories
DNA Fragmentation	/	✓ .	USB 3.0 CCD Camera
Leukocytes	1	✓ •	23" touch screen all-in-one computer with stylus
Antisperm - Antibody		NA •	barcode reader & printer
α - Glycosidase	NA	✓ .	10" tablet - windows OS
Fructose	NA	✓ .	colour inkjet printer
Zinc	NA	✓ .	androchem analyser - optiona (with R32 cable)



a non-spermicidal semen collection pouch made of biologically inert polyurethane

maintains sperm motility & viability






SpermTrans is an Intrauterine Insemination Catheter made up of medical grade Polyolefin


working designs suits uterine anatomy



Product Code SP/PAD/001







Andro Table

ergonomically designed table which makes lowest foot print andrology lab



to further add to our 360° solutions

We also provide

Hands-on Training

on Basic & Advance Lab Andrology Techniques

thus we are with you 360 degrees to make your lab complete



www.spermprocessor.com

Sperm Processor Pvt. Ltd.

6, Welcomenagar, Garkheda, Aurangabad - 431005, India

Ph.: +91 240 6603800 Fax: +91 240 2341694

Email: info@spermprocessor.com